


How to use and washing instructions for
No. 98W Easy Grip Pill Crusher


HOW TO USE




Put pill into crusher and turn

Twist on/off lid for storage


Washing by hand



Open Pill Crusher and rinse under tap



Carefully dry with a tissue



Please issue to users

W+W MEDSYSTEMS
LIMITED

2-4 Lockwood House, Lockwood Park
Brewery Drive, Huddersfield HD4 6EN

Tel: 01484 667822 Fax: 01484 667817


Email: info@wwmed.co.uk

Website: www.wwmed.co.uk

Sterilisation

Steam time and temperature

18 minutes at 121°C or 3-4 minutes at 134°C



Contains small parts that could be a choking hazard for children or any individuals who have a tendency to place objects in their mouths

This device complies with EC 1207/2001 and is Year 2000 and Euro compliant